

# Clostridium K-Set


**Rapid diagnostic test for *in vitro* detection of *Clostridium difficile* antigen in stool specimens**

- Safe
- Easy and ready to use
- Complete set of sampling devices available
- Accurate
- Easy to interpret



[www.corisbio.com](http://www.corisbio.com)

# Clostridium K-SeT



## Performance compared to Culture (n=318)

|             |       | (95% confidence interval) |
|-------------|-------|---------------------------|
| Specificity | 93%   | (88.8 to 95.7%)           |
| Sensitivity | 100%  | (93.9 to 100%)            |
| PPV         | 81.5% | (71.8 to 88.6%)           |
| NPV         | 100%  | (97.9 to 100%)            |
| Accuracy    | 94.7% |                           |

## For more technical information

[client.care@corisbio.com](mailto:client.care@corisbio.com)

## Ordering Information

| Product name      | Description                              | Storage     | Code   |
|-------------------|--|-------------|--------|
| Clostridium K-SeT | 20 cassettes, buffer                     | 4°C to 30°C | K-1220 |
| Clostridium K-SeT | 20 cassettes, 20 faecal sampling systems | 4°C to 30°C | K-1520 |


### All diagnostics available in Coris BioConcept catalogue:

Adenovirus, Adenovirus (Respiratory), Adenovirus 40 and 41, *Clostridium difficile*, Bovine Coronavirus, *Cryptosporidium*, *Escherichia coli* O157, *Escherichia coli* K99, *Giardia*, *Helicobacter pylori*, Influenza A and B viruses, *Legionella pneumophila*, Respiratory Syncytial Virus, Rotavirus.

Mefloquine, Proguanil.

*Leishmania spp.*, *Trypanosoma cruzi*, *P. aeruginosa* efflux.

Distributed by:



Science Park CREALYS  
Rue Jean Sonet 4A - 5032 Gembloux - Belgium  
Ph.: +32 (0) 81 719 910 - Fx.: +32 (0)81 719 919  
[sales@corisbio.com](mailto:sales@corisbio.com) - [www.corisbio.com](http://www.corisbio.com)

MADE IN BELGIUM

Fo 7.5-34-20/01  
January 2012