

0157 Coli-Strip


Rapid diagnostic test for *in vitro* detection of *Escherichia coli* O157 bacteria in stool specimens (after broth enrichment)

- Easy to use
- Easy to read
- Specific
- Sensitive
- Reliable
- Rapid (15')



www.corisbio.com

O157 Coli-Strip



Performance

Limit of detection

The O157 Coli-Strip's limit of detection for *E. coli* O157 is 10^4 CFU/mL.

Interference

No cross-reactivity with enteric bacteria, parasites and/or viruses.

For more technical information

client.care@corisbio.com


Ordering Information

Product Name	Description	Storage	Code
O157 Coli-Strip	20 tests, buffer	4°C to 30°C	C-1011
O157 Coli Uni-Strip	10 individual strips, buffer sample taking material (test tubes, loops, tube holder)	4°C to 30°C	K-1511

All diagnostics available in Coris BioConcept catalogue :

Adenovirus, Adenovirus (Respiratory), Adenovirus 40 and 41, *Clostridium difficile*, Bovine Coronavirus, *Cryptosporidium*, *Escherichia coli* O157, *Escherichia coli* K99, Giardia, *Helicobacter pylori*, Influenza A and B viruses, *Legionella pneumophila*, Respiratory Syncytial Virus, Rotavirus, Mefloquine, Proguanil, *Leishmania spp.*, *Trypanosoma cruzi*, *P. aeruginosa* efflux.

Distributed by:



Science Park CREALYS
Rue Jean Sonet 4A - 5032 Gembloux - Belgium
Ph.: +32 (0) 81 719 910 - Fx.: +32 (0)81 719 919
sales@corisbio.com - www.corisbio.com

MADE IN BELGIUM

Fo 7.5-32-11/04
January 2012