


RSV K-Set

Rapid diagnostic test for *in vitro* detection of Respiratory Syncytial Virus in nasopharyngeal specimens

- Safe
- Easy and ready to use
- Complete set of sampling devices available
- Accurate
- Easy to interpret



www.corisbio.com



Performance compared to

	IF (n=100)	DFA (n=236)
Specificity	100% (50/50)	98.3% (119/121)
Sensitivity	96% (48/50)	92.2% (106/115)
PPV	100% (48/48)	98.1% (106/108)
NPV	96.1% (50/52)	93% (119/128)
Accuracy	98% (98/100)	95.3% (225/236)

For more technical information

client.care@corisbio.com


Ordering Information

Product Name	Description	Storage	Code
RSV K-SeT	20 cassettes, buffer	4°C to 30°C	K-1206
RSV K-SeT, with collection set	20 cassettes, buffer, 20 swabs	4°C to 30°C	K-1506
RSV Respi-Strip	25 tests, buffer	4°C to 30°C	C-1006

All diagnostics available in Coris BioConcept catalogue :

Adenovirus, Adenovirus (Respiratory), Adenovirus 40 and 41, *Clostridium difficile*, Bovine Coronavirus, *Cryptosporidium*, *Escherichia coli* O157, *Escherichia coli* K99, Giardia, *Helicobacter pylori*, Influenza A and B viruses, *Legionella pneumophila*, Respiratory Syncytial Virus, Rotavirus, Mefloquine, Proguanil, *Leishmania spp.*, *Trypanosoma cruzi*, *P. aeruginosa* efflux.

Distributed by:



Science Park CREALYS
Rue Jean Sonet 4A - 5032 Gembloux - Belgium
Ph.: +32 (0) 81 719 910 - Fx.: +32 (0)81 719 919
sales@corisbio.com - www.corisbio.com

MADE IN BELGIUM